



EONICAAESTHETICS
ETERNAL APPEAL

Artos


PICOSECOND - ERASE. RESTORE. REJUVENATE

ARTOS BEAUTY THROUGH INNOVATION



EONICAESTHETICS ETERNAL APPEAL

MECHANISM OF ACTION

The powerful Active Q-Switch YAG laser, featuring an ultra-short pulse width of 2 nanoseconds, is an exceptionally effective technique for eliminating both natural and artificial pigmentation. The laser's energy is absorbed by the skin's chromophores, leading to a photoacoustic effect that breaks down the pigment while preserving the surrounding tissues. This process ultimately results in the gradual fading and clearing of pigmentation as the body's immune system responds.




Advantages




01/ Five treatment modes

There are five available treatment modes, each ensuring the ideal pulse width and energy output for various indications, resulting in optimal outcomes and minimal recovery time!!




02/ Built-in clinical database

The integrated clinical database offers recommended parameters that are frequently used for easy reference, simplifying the process for beginners to get started!




03/ Four independent control systems with integrated power supply.

The combination of an integrated power supply and four capacitors guarantees not only the highest stability in power output but also serves as a crucial distinguishing feature of the genuine Pico Qswitch yag Laser (most competitors utilize just one capacitor)!



04/ Simple screw-off external bottle filters

The convenient, removable filter bottle purifies distilled water, creating the ideal environment for advanced precision laser sources, which significantly extends their lifespan!



05/ Ultimate energy yield

When using LONG PULSE, energy is at 2500mj;
Energy is 800mj in picoseconds.



06/ Interlock

Introducing interlock system, providing an additional security layer for clients.

ARTOS BEAUTY THROUGH INNOVATION



Tattoos removal



Freckle removal



Hair removal



Acne clearance



Improve varicose
veins



Facial rejuvenation



Onychomycosis



Brief pulse duration, Intense immediate energy

Ultra-short pulse widths deliver a more intense impact on the pigment, breaking it down into smaller particles to enhance metabolism and eliminate it from the body.



Diverse wavelengths

Thoroughly tackles all deep, stubborn pigments and harmless pigmented lesions while addressing various skin concerns.




Precise

Innovative technology ensures 100% energy precision to control both the depth and intensity of treatment.

EA





- Customized **Korean Seven-Joint Arm**
- **High barrel light guide arm:** rate of energy loss $\leq 8\%$
- Modern design
- **360°** of rotation, ensures constant energy laser power
- It is more **shock-resistant** than existing products

APPLICATIONS

PICOSECOND

It can be applied to most of the spots, birthmarks or tattoos, skin rejuvenation treatment, tattoo!

PTP

Melasma, Moles, Tattoo, Dark Eye Circles, Pores, Hyperpigmentation

MULTI PULSE


Mainly used for chloasma, Stubborn stain and pigments

LONG PULSE

Vascular, Scars, Acne, Rosacea/dif fuse redness, Onychomycosis, Rejuvenation, firming, whitening, shrinking pores, Remove fine wrinkles, Heating dermal cells, hair removal, improve Varicose veins

LONG PULSE Pro

Vascular, Scars, Acne, Rosacea/dif fuse redness, Onychomycosis, Rejuvenation, firming, whitening, shrinking pores, Remove fine wrinkles, Heating dermal cells, hair removal, improve Varicose veins



Acne Scars




depressed scarring




dull skin



uneven skin tone




lack of collagen



stubborn tattoos



Focusing treatment head

wavelength: 1064nm, 532nm




Honeycomb treatment head

wavelength: 755nm



Before and After



Artos PICOSECOND - ERASE. RESTORE. REJUVENATE



TECHNICAL SPECIFICATIONS

Beam profile	Flat-top mode				
Wavelength	1064nm,532nm,755nm				
Treatment modes	Picosecond / PTP / Multi pulse / Long pulse / Long pulse Pro				
Pulse energy	1064NM	Model	Minimum	Maximum	Stepping
		Picosecond	100mj	800mj	50mj
		PTP	100mj	1200mj	50mj
		Multi Pulse	100mj	1500mj	50mj
		Long Pulse	200mj	2500mj	100mj
		Long Pulse Pro	200mj	5000mj	100mj
	532NM	Picosecond	50mj	400mj	25mj
		PTP	50mj	600mj	25mj
		Multi Pulse	50mj	750mj	25mj
		Long Pulse	/	/	/
		Long Pulse Pro	/	/	/
Spot size	2-10mm Adjustable				
Aiming beam	Diode 650nm<5mw (Red)				
Repetition rate	1-10HZ adjustable				
Package dimension and weight	100*53*115cm (101KG) 21*44*106 (13KG)				
Power supply	AC110V/230V, 50/60HZ				

EONICAAESTHETICS
ETERNAL APPEAL